

Product Datasheet

Mouse OX40 / TNFRSF4 / CD134 Protein (orb257730)

Description Mouse OX40, His Tag is expressed from human 293 cells (HEK293). It contains AA Val 20 - Pro 211 (Acc...

Reactivity Mouse

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target OX40

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag C-10×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 23.2 kDa. The protein migrates as 40-55 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_035789.1

Purity 90%


Source Mouse OX40, His Tag (orb257730) is expressed from human 293 cells (HEK293). It contains AA Val 20 - Pro 211 (Accession # P47741-1).

MW 23.2 kDa


Uniprot ID [P47741](#)

NCBI [NP_035789.1](#)

Expiration Date 6 months from date of receipt.



SDS-PAGE analysis of Mouse OX40 protein



Line graph illustrates about the sensitivity...