


## Product Datasheet

### Human CD44 Protein (orb257773)

<b>Description</b>	Human CD44 Protein, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Gln 21 - Pro 2...
<b>Reactivity</b>	Human
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	CD44
<b>Preservatives</b>	Tris with Glycine, Arginine and NaCl, pH7.5
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C
<b>Tag</b>	C-Fc
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 48.7 kDa. The protein migrates as 55-70 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Protein Sequence</b>	AAH04372
<b>Purity</b>	95%
<b>Source</b>	Human CD44 Protein, Fc Tag (orb257773) is expressed from human 293 cells (HEK293). It contains AA Gln 21 - Pro 220 (Accession # AAH04372).
<b>MW</b>	48.7 kDa
<b>Uniprot ID</b>	<b>P16070</b>
<b>NCBI</b>	<b>AAH04372</b>
<b>Expiration Date</b>	6 months from date of receipt.



SDS-PAGE analysis of Human CD44 protein