


Product Datasheet

Human PLGF / PGF Protein (orb257774)

Description	Human PIGF-2 Protein, His Tag is expressed from human 293 cells (HEK293). It contains AA Leu 19 - Ar...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	PLGF
Preservatives	100 mM Acetic Acid, pH3.0
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 18.2 kDa. The protein migrates as 28-33 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_002623.2
Purity	95%
Source	Human PLGF, His Tag (orb257774) is expressed from human 293 cells (HEK293). It contains AA Leu 19 - Arg 170 (Accession # NP_002623.2).
MW	18.2 kDa
Uniprot ID	P49763
NCBI	NP_002623.2
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human PLGF protein