


Product Datasheet

Human plgR Protein (orb257776)

Description	Human plgR, His Tag is expressed from human 293 cells (HEK293). It contains AA Lys 19 - Arg 638 (Acc...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	plgR
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-6×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 68.7 kDa. The protein migrates as 90-100 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAI10495.1
Purity	95%
Source	Human plgR, His Tag (orb257776) is expressed from human 293 cells (HEK293). It contains AA Lys 19 - Arg 638 (Accession # AAI10495.1).
MW	68.7 kDa
Uniprot ID	P01833
NCBI	AAI10495.1
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human plgR protein