


## Product Datasheet

### Mouse ROR1 Protein (orb257798)

<b>Description</b>	Mouse ROR1, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Gln 30 - Glu 403 (Acce...
<b>Reactivity</b>	Mouse
<b>Endotoxins</b>	1.0 EU per µg
<b>Conjugation</b>	Unconjugated
<b>Target</b>	ROR1
<b>Preservatives</b>	50 mM Tris, 100 mM Glycine, 25 mM Arginine, 150 mM NaCl, pH7.5
<b>Form/Appearance</b>	Powder
<b>Storage</b>	-20°C
<b>Tag</b>	C-Fc
<b>Note</b>	For research use only
<b>Application notes</b>	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 68.4 kDa. The protein migrates as 76-96 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
<b>Protein Sequence</b>	Gln 30 - Glu 403
<b>Purity</b>	95%
<b>Source</b>	Mouse ROR1, Fc Tag (orb257798) is expressed from human 293 cells (HEK293). It contains AA Gln 30 - Glu 403 (Accession # Q9Z139).
<b>MW</b>	68.4 kDa
<b>Uniprot ID</b>	<a href="#">Q9Z139</a>
<b>Expiration Date</b>	6 months from date of receipt.



SDS-PAGE analysis of Mouse ROR1 protein