



Product Datasheet

Human ROR2 / NTRKR2 Protein (orb257799)

Description	Human ROR2, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Glu 34 - Gly 403 (Acce...
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	ROR2
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 68 kDa. The protein migrates as 70-90 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAI30523.1
Purity	95%
Source	Human ROR2, Fc Tag (orb257799) is expressed from human 293 cells (HEK293). It contains AA Glu 34 - Gly 403 (Accession # A1L4F5-1).
MW	68 kDa
Uniprot ID	A1L4F5
NCBI	AAI30523.1
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human ROR2 protein



Line graph illustrates about the Ag-Ab r...