

## Product Datasheet

### Human R-Spondin 1 / RSP01 Protein (orb257807)

**Description** Human R-Spondin 1, His Tag is expressed from human 293 cells (HEK293). It contains AA Ser 21 - Ala 2...

**Reactivity** Human

**Endotoxins** 1.0 EU per µg

**Conjugation** Unconjugated

**Target** R-Spondin 1

**Preservatives** PBS, pH7.4

**Form/Appearance** Powder

**Storage** -20°C

**Tag** C-10×His

**Note** For research use only

**Application notes** This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 28.7 kDa. The protein migrates as 38-42 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

**Protein Sequence** NP\_001033722.1

**Purity** 95%


**Source** Human R-Spondin 1, His Tag (orb257807) is expressed from human 293 cells (HEK293). It contains AA Ser 21 - Ala 263 (Accession # Q2MKA7-1).

**MW** 28.7 kDa

**Uniprot ID** [Q2MKA7](#)

**NCBI** [NP\\_001033722.1](#)

**Expiration Date** 6 months from date of receipt.



SDS-PAGE  
analysis of  
Human R-  
Spondin 1  
p...