

Product Datasheet

Human TGF-beta RII / TGFBR2 Protein (orb257866)

Description Human TGFBR2, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Thr 23 - Asp 159 (Ac...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target TGF-beta RII

Preservatives Tris with Glycine, Arginine and NaCl, pH7.5

Form/Appearance Powder

Storage -20°C

Tag C-Fc

Note For research use only

Application notes This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 42.1 kDa. The protein migrates as 55-66 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_003233

Purity 95%


Source Human TGFBR2, Fc Tag (orb257866) is expressed from human 293 cells (HEK293). It contains AA Thr 23 - Asp 159 (Accession # NP_003233).

MW 42.1 kDa

Uniprot ID **P37173**

NCBI **NP_003233**

Expiration Date 6 months from date of receipt.



SDS-PAGE
analysis of
Human
TGF-beta
RII ...