

Product Datasheet

Human TNFR1 / CD120a / TNFRSF1A Protein (orb257876)

Description

Human TNFR1, His Tag is expressed from human


Reactivity

(Ac...

Human

Endotoxins

1.0 EU per µg

Conjugation

Unconjugated

Target

TNFR1

Preservatives

PBS, pH7.4

Form/Appearance

Powder

Storage

-20°C

Tag

C-8×His

Note

For research use only

Application notes

This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 22.3 kDa. The protein migrates as 28-40 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to different glycosylation.

Protein Sequence

NP_001056

Purity

90%

Source

Human TNFR1, His Tag (orb257876) is expressed from human 293 cells (HEK293). It contains AA Ile 22 - Thr 211 (Accession # NP_001056).

MW

22.3 kDa

Uniprot ID
P19438
NCBI
NP_001056
Expiration Date

6 months from date of receipt.

SDS-PAGE
analysis of
Human
TNFR1
protein