

Product Datasheet

Human TRAIL R1 / DR4 / TNFRSF10A Protein (orb257894)

Description

Human TRAIL R1, His Tag is expressed from



Reactivity Asn 239 ...

Reactivity Human

Endotoxins 1.0 EU per µg

Conjugation Unconjugated

Target TRAIL R1

Preservatives PBS, pH7.4

Form/Appearance Powder

Storage -20°C

Tag C-10×His

Note For research use only

Application notes This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 24.8 kDa. The protein migrates as 20-23 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.

Protein Sequence NP_003835.3

Purity 90%


Source Human TRAIL R1, His Tag (orb257894) is expressed from human 293 cells (HEK293). It contains AA Ala 24 - Asn 239 (Accession # O00220-1).

MW 24.8 kDa

Uniprot ID [O00220](https://www.uniprot.org/entry/O00220)

NCBI [NP_003835.3](https://www.ncbi.nlm.nih.gov/nuccore/NP_003835.3)

Expiration Date 6 months from date of receipt.



SDS-PAGE analysis of Human TRAIL R1 prot...