


Product Datasheet

Human Vasorin / VASN Protein (orb257915)

Description	Human Vasorin, His Tag is expressed from human 293 cells (HEK293). It contains AA Cys 24 - Pro 575 (...)
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Vasorin
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 61.6 kDa. The protein migrates as 80-100 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAH68575
Purity	95%
Source	Human Vasorin, His Tag (orb257915) is expressed from human 293 cells (HEK293). It contains AA Cys 24 - Pro 575 (Accession # AAH68575).
MW	61.6 kDa
Uniprot ID	Q6EMK4
NCBI	AAH68575
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Human Vasorin prote...