

Product Datasheet

**Ebolavirus EBOV (subtype Zaire, strain Kikwit-95)
Envelope Glycoprotein 1 (GP1) Protein
(orb257933)**

Description	Ebolavirus EBOV (subtype Zaire, strain from human 293 cells (HE...))	
Reactivity	Virus	
Endotoxins	1.0 EU per µg	SDS-PAGE analysis of Viral Glycoprotein ...
Conjugation	Unconjugated	
Target	Glycoprotein / GP (virus)	
Preservatives	PBS, pH7.4	
Form/Appearance	Powder	
Storage	-20°C	
Tag	C-8×His	
Note	For research use only	
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 52.0 kDa. The protein migrates as 95-120 kDa under reducing (R) condition (SDS-PAGE).	
Protein Sequence	AAQ55048.1	
Purity	95%	
Source	Ebolavirus EBOV (subtype Zaire, strain Kikwit-95) GP1, His Tag (orb257933) is expressed from human 293 cells (HEK293). It contains AA Ile 33 - Arg 501 (Accession # AAQ55048.1).	
MW	52.0 kDa	
Uniprot ID	P87666	
NCBI	AAQ55048.1	
Expiration Date	6 months from date of receipt.	