


Product Datasheet

Rhesus macaque CD40 / TNFRSF5 Protein (orb257939)

Description	Rhesus macaque CD40, His Tag is expressed from human 293 cells (HEK293). It contains AA Glu 21 - Arg...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	CD40
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 21.1 kDa. The protein migrates as 30-34 kDa under reducing (R) condition (SDS-PAGE).
Protein Sequence	NP_001252791.1
Purity	95%
Source	Rhesus macaque CD40, His Tag (orb257939) is expressed from human 293 cells (HEK293). It contains AA Glu 21 - Arg 193 (Accession # NP_001252791.1).
MW	21.1 kDa
Uniprot ID	NP_001252791.1
NCBI	NP_001252791.1
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Monkey CD40 protein