


Product Datasheet

Rhesus macaque EGF R Protein (orb257945)

Description	Rhesus macaque EGF R, Fc Tag is expressed from human 293 cells (HEK293). It contains AA Leu 25 - Ser...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	EGF R
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 95.3 kDa. The protein migrates as 116-120 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Leu 25 - Ser 645
Purity	95%
Source	Rhesus macaque EGF R, Fc Tag (orb257945) is expressed from human 293 cells (HEK293). It contains AA Leu 25 - Ser 645 (Accession # P55245).
MW	95.3 kDa
Uniprot ID	P55245
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Monkey EGF R protei...