


Product Datasheet

Rhesus macaque EGF R Protein (orb257946)

Description	Rhesus macaque EGF R, His Tag is expressed from human 293 cells (HEK293). It contains AA Leu 25 - Ser...
Reactivity	Monkey
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	EGF R
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 70.5 kDa. The protein migrates as 80-100 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	Leu 25 - Ser 645
Purity	95%
Source	Rhesus macaque EGF R, His Tag (orb257946) is expressed from human 293 cells (HEK293). It contains AA Leu 25 - Ser 645 (Accession # P55245).
MW	70.5 kDa
Uniprot ID	P55245
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Monkey EGF R protei...