



Product Datasheet

Mouse LDL R Protein (orb257948)

Description	Mouse LDL R, His Tag is expressed from human 293 cells (HEK293). It contains AA Ala 22 - Arg 790 (Ac...
Reactivity	Mouse
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	LDL R
Preservatives	PBS, pH7.4
Form/Appearance	Powder
Storage	-20°C
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 86.2 kDa. The protein migrates as 90-130 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	NP_034830
Purity	95%
Source	Mouse LDL R, His Tag (orb257948) is expressed from human 293 cells (HEK293). It contains AA Ala 22 - Arg 790 (Accession # NP_034830).
MW	86.2 kDa
Uniprot ID	O88307
NCBI	NP_034830
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Mouse LDL R protein



Line graph illustrates about the sensitivity...