


Product Datasheet

HCMV Glycoprotein B (gB) Protein (orb257973)

Description	HCMV Glycoprotein B (gB) Protein, Fc Tag is expressed from human 293 cells (HEK293). It contains AA ...
Reactivity	HCMV
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	Glycoprotein B / gB
Preservatives	Tris with Glycine, Arginine and NaCl, pH7.5
Form/Appearance	Powder
Storage	-20°C
Tag	C-Fc
Note	For research use only
Application notes	This protein carries a human IgG1 Fc tag at the C-terminus. The protein has a calculated MW of 104.0 kDa. The protein migrates as 130-180 kDa under reducing (R) condition (SDS-PAGE) due to glycosylation.
Protein Sequence	AAA45920.1
Purity	90%
Source	HCMV Glycoprotein B (gB) Protein, Fc Tag (orb257973) Val 23 - Lys 700 (Accession # P13201-1) with furin cleavage site mutated from 'RTKR' to 'TTQT', was produced in human 293 cells (HEK293) at ACROBiosystems.
MW	104.0 kDa
Uniprot ID	P13201
NCBI	AAA45920.1
Expiration Date	6 months from date of receipt.



SDS-PAGE analysis of Viral Glycoprotein ...