



Product Datasheet

SYNJ2BP antibody (orb25864)


Description	SYNJ2BP antibody
Species/Host	Rabbit
Reactivity	Human, Rat
Conjugation	Unconjugated
Tested Applications	ELISA, IHC, IP, WB
Immunogen	Recombinant Human Synaptojanin-2-binding protein (1-117AA)
Target	SYNJ2BP
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P57105
Dilution Range	WB:1:500-5000, IHC-P:1:20-200, IP:1:200-2000
Expiration Date	12 months from date of receipt.



Western blot analysis of MCF7(lane 1), U...



Immunohistochemical staining of human te...



Immunoprecipitation analysis of Rat kidn...