

Product Datasheet

Anti-TRAP/CD40L/Cd40Ig Antibody (orb2599415)

| | |
|----------------------------|--|
| Description | Anti-TRAP/CD40L/Cd40Ig Antibody. Tested in ELISA, IHC, WB applications. This antibody reacts with Mouse, Rat. |
| Species/Host | Rabbit |
| Reactivity | Mouse, Rat |
| Conjugation | PE |
| Tested Applications | ELISA, IHC, WB |
| Immunogen | E.coli-derived rat TRAP/CD40L/Cd40Ig recombinant protein (Position: M112-L260). |
| Form/Appearance | Lyophilized |
| Concentration | Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml. |
| Storage | At -20°C for one year from date of receipt. After reconstitution, at 4°C for one month. It can also be aliquotted and stored frozen at -20°C for six months. Avoid repeated freezing and thawing. |
| Note | For research use only |
| Application notes | Western blot, 0.25-0.5 µg/ml, Mouse, Rat Immunohistochemistry(Paraffin-embedded Section), 2-5 µg/ml, Rat ELISA, 0.1-0.5 µg/ml, Rat. Adding 0.2 ml of distilled water will yield a concentration of 500 µg/ml |
| Isotype | Rabbit IgG |
| Clonality | Polyclonal |
| MW | 29 kDa |
| Uniprot ID | Q9Z2V2 |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)