

Product Datasheet

Donkey IgG (H&L) - Cy5.5 (orb2624466)

Description	Donkey IgG
Species/Host	Rabbit
Reactivity	Donkey
Conjugation	Cy5.5
Tested Applications	ELISA, FC, ICC, IHC
Target	Donkey IgG (H&L) - Cy5.5
Concentration	batch dependent
Storage	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
Note	For research use only
Isotype	Rabbit
Clonality	Polyclonal
Dilution Range	IHC: 1/100 - 1/1000; ICC: 1/100 - 1/1000; FCM: 1/1000 - 1/4000; ELISA: Use at an assay dependent concentration
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)