

## Product Datasheet

### Donkey IgG (H&L) - Cy5 (orb2624495)

<b>Description</b>	Donkey IgG
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Donkey
<b>Conjugation</b>	Cy5
<b>Tested Applications</b>	ELISA, FC, ICC, IHC
<b>Target</b>	Donkey IgG (H&L) - Cy5
<b>Concentration</b>	batch dependent
<b>Storage</b>	Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit
<b>Clonality</b>	Polyclonal
<b>Dilution Range</b>	IHC: 1/100 - 1/1000; ICC: 1/100 - 1/1000; FCM: 1/1000 - 1/4000; ELISA: Use at an assay dependent concentration
<b>Expiration Date</b>	12 months from date of receipt.

---

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)