

Product Datasheet

Human IgG Fab - Texas Red (orb2624572)

| | |
|----------------------------|--|
| Description | Human IgG Fab |
| Species/Host | Mouse |
| Reactivity | Human |
| Conjugation | Texas Red |
| Tested Applications | ELISA, FC, ICC, IHC |
| Target | Human IgG Fab - Texas Red |
| Preservatives | 1 mg/ml. Liquid in 0.01M Phosphate Buffered Saline, pH 7.2, containing 1% rAlbumin, 50% glycerol, 0.02% Sodium Azide |
| Concentration | 1mg/ml |
| Storage | Shipped at 4°C. Upon delivery aliquot and store at -20°C for one year. Avoid freeze/thaw cycles. |
| Note | For research use only |
| Isotype | Mouse |
| Clonality | Monoclonal |
| Dilution Range | IHC: 1/100 - 1/1000; ICC: 1/100 - 1/1000; FCM: 1/1000 - 1/4000; ELISA: Use at an assay dependent concentration |
| Expiration Date | 12 months from date of receipt. |

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)