

Product Datasheet

Tyrosine Hydroxylase (Ab-19) Antibody (orb2627803)

Description	Synthesized non-phosphopeptide derived from human Tyrosine Hydroxylase around the phosphorylation site of serine 19 (A-V-S(p)-E-Q).
Species/Host	Rabbit
Reactivity	Human, Mouse, Rat
Conjugation	Unconjugated
Tested Applications	IF, IHC, WB
Immunogen	Synthesized non-phosphopeptide derived from human Tyrosine Hydroxylase around the phosphorylation site of serine 19 (A-V-S(p)-E-Q).
Target	Tyrosine Hydroxylase (Ab-19)
Preservatives	Purified antibody in PBS with 0.05% sodium azide.
Concentration	Approximately 1mg/ml (varies from batch to batch). Please inquire for precise concentration.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG
Clonality	Polyclonal
MW	56kDa
Uniprot ID	P07101
Entrez	7054

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Dilution Range WB: 1/500-1/3000; IHC: 1/50-1/100; ICC: 1/100-1/500

Expiration Date 12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)