

## Product Datasheet

### MDM2(Ab-166) Antibody (orb2627866)

<b>Description</b>	Peptide sequence around aa.164~168 (A-T-S-E-T) derived from Human MDM2.
<b>Species/Host</b>	Rabbit
<b>Reactivity</b>	Human, Monkey, Mouse, Rat
<b>Conjugation</b>	Unconjugated
<b>Tested Applications</b>	ICC, IHC, WB
<b>Target</b>	MDM2(Ab-166)
<b>Preservatives</b>	Supplied at 1.0mg/mL in phosphate buffered saline (without Mg <sup>2+</sup> and Ca <sup>2+</sup> ), pH 7.4, 150mM NaCl, 0.02% sodium azide and 50% glycerol.
<b>Concentration</b>	1.0mg/mL
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Isotype</b>	Rabbit IgG
<b>Clonality</b>	Polyclonal
<b>MW</b>	90kDa
<b>Uniprot ID</b>	<a href="#">Q00987</a>
<b>Entrez</b>	<a href="#">4193</a>
<b>Dilution Range</b>	WB: 1/500-1/1000 ; IHC: 1/50-1/100 ; ICC: 1/100-1/200
<b>Expiration Date</b>	12 months from date of receipt.