

Product Datasheet

Eribulin Rabbit Recombinant Monoclonal Antibody (orb2645274)

Description	Monoclonal Anti-Eribulin Antibody, Rabbit IgG (1M2B11) is a Rabbit monoclonal antibody recombinantly expressed from HEK293 cells.
Species/Host	Rabbit
Reactivity	Chemical
Conjugation	Unconjugated
Tested Applications	ELISA
Target	Eribulin
Preservatives	PBS, pH 7.4
Form/Appearance	Lyophilized from 0.22 µm filtered solution in PBS, pH 7.4 with trehalose as protectant. Contact us for customized product form or formulation.
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only
Isotype	Rabbit IgG Rabbit Kappa
Clonality	Recombinant
Clone Number	1M2B11
Antibody Type	Recombinant Antibody
Purity	> 95% as determined by SDS-PAGE. Protein A purified / Protein G purified.
Dilution Range	ELISA: 0.2-125 ng/mL

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)



Explore. Bioreagents.

Biorbyt.com

Expiration Date

12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)