

## Product Datasheet

### Direct qPCR ProbesMaster highROX (orb2646469)

<b>Description</b>	Robust real-time PCR master mix with highROX for amplification directly from blood, swabs or tissue
<b>Storage</b>	store at -20 °C. avoid freeze/thaw cycles,. stable at 4 °C for up to 4 weeks
<b>Note</b>	For research use only
<b>Expiration Date</b>	12 months from date of receipt.

---

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)