

Product Datasheet

AffiniPure Mouse Anti-Rabbit IgG (H+L) secondary antibody, FITC Conjugated (orb2646591)

Description	AffiniPure Mouse Anti-Rabbit IgG (H+L) secondary antibody, FITC Conjugated
Species/Host	Mouse
Reactivity	Rabbit
Conjugation	FITC
Tested Applications	FC, IF
Immunogen	Rabbit IgG (whole molecular)
Form/Appearance	Concentrated, Liquid
Concentration	1 mg/ml
Storage	At -20°C for one year from date of receipt. Avoid repeated freezing and thawing. Protect from light.
Note	For research use only
Application notes	Immunofluorescence, 10-20 µg/mL Flow Cytometry, 10-20 µg/mL
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)