

Product Datasheet

PE Anti-Mouse CD40(FGK4.5/FGK45) (orb2648339)

Description	PE Anti-Mouse CD40(FGK4.5/FGK45)
Species/Host	Rat
Reactivity	Mouse
Conjugation	PE
Tested Applications	FC
Form/Appearance	liquid
Storage	PBS with 0.1% sodium azide, 1%BSA, pH 7.3, 2-8°C for 12 months (Avoid repeated freeze / thaw cycles.)
Note	For research use only
Isotype	IgG2a, kappa
Clonality	Monoclonal
Clone Number	FGK4.5/FGK45
Uniprot ID	P27512
Dilution Range	Volume per test: 5µL. Each lot of this antibody is quality control tested by flow cytometric analysis. The amount of the reagent is suggested to be used 5 µL of antibody per test (million cells in 100 µL staining volume or per 100 µL of whole blood). Please check your vial before the experiment. Since applications vary, the appropriate dilutions must be determined for individual use.
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)