

Product Datasheet

Streptavidin Conjugated with Fluorescein (FITC) (orb2652868)

Description	Streptavidin is a 55 kDa (subunit MW 14 kDa) biotin-binding protein isolated from <i>Streptomyces avidini</i> . Streptavidin is superior to avidin, because it does not contain carbohydrate like avidin and has no net charge at neutral pH. Streptavidin~biotin system is routinely used in Immunohistochemistry (IHC).
Tested Applications	IHC
Storage	2-8 °C, in dark
Note	For research use only
Application notes	To be used with Biotin conjugated proteins.
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)