

Product Datasheet

Senkyunolide N (orb2653517)

Description	Senkyunolide N exhibits prostaglandin H endoperoxide synthase-I (COX-I) and endoperoxide synthase-II (COX-II) Inhibitory activity at pH 7 and at a concentration of 250 µg ml ⁻¹ .
Storage	Storage temperature: -20°C. Stability: ≥ 2 years
Note	For research use only
Formula	C ₁₂ H ₁₈ O ₄
Purity	> 98% (HPLC)
MW	226.3
CAS Number	140694-58-2
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)