

Product Datasheet

Polygalasaponin XXXI (orb2654545)

Description	Onjisaponin F(Polygalasaponin XXXI) acts non-competitively against cyclic AMP phosphodiesterase, like papaverine; it also exhibits a prolongation effect on hexobarbital sleeping time in mice.
Target	cAMP
Storage	Storage temperature: -20°C. Stability: ≥ 2 years
Note	For research use only
Formula	$C_{75}H_{112}O_{36}$
Purity	> 98% (HPLC)
MW	1589.7
CAS Number	79103-90-5
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)