

Product Datasheet

Isoferulic acid (orb2655027)

Description	Isoferulic acid is an effective natural antioxidant in both lipid and aqueous media, it may be a new promising anti-glycation agent for the prevention of diabetic complications via inhibition of advanced glycation end products (AGEs) formation and oxidation-dependent protein damage.
Storage	Storage temperature: -20°C. Stability: ≥ 2 years
Note	For research use only
Formula	$C_{10}H_{10}O_4$
Purity	> 98% (HPLC)
MW	194.2
CAS Number	25522-33-2
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+441223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)