

Product Datasheet

SPR741 (orb2655224)

Description	SPR741, formerly NAB741 (12), is a cationic peptide derived from polymyxin B and is one such potentiator molecule that is under development for the treatment of serious Gram-negative bacterial infections.
Storage	Storage temperature: -20°C. Stability: ≥ 2 years
Note	For research use only
Formula	$C_{44}H_{73}N_{13}O_{13}$
Purity	> 98% (HPLC)
MW	992.1
CAS Number	1179330-52-9
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)