

Product Datasheet

(R)-Tegoprazan (orb2655256)

Description	(R)-Tegoprazan (example 3), a benzimidazole derivative, is a potent kidney H ⁺ /K ⁺ -ATPase inhibitor with an IC ₅₀ of 98 nM of canine kidney Na ⁺ /K ⁺ -ATPase.
Storage	Storage temperature: -20°C. Stability: ≥ 2 years
Note	For research use only
Formula	C ₂₀ H ₁₉ F ₂ N ₃ O ₃
Purity	> 98% (HPLC)
MW	387.4
CAS Number	942195-56-4
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)