

Product Datasheet

SR-4835 (orb2655432)

Description	SR-4835 is a highly selective dual inhibitor of CDK12 and CDK13(CDK12: IC50=99 nM Kd=98 nM; CDK13: Kd=4.9 nM) which disables triple-negative breast cancer (TNBC) cells.
Target	CDK
Storage	Storage temperature: -20°C. Stability: ≥ 2 years
Note	For research use only
Formula	$C_{21}H_{20}Cl_2N_{10}O$
Purity	> 98% (HPLC)
MW	499.4
CAS Number	2387704-62-1
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)