

Product Datasheet

Fingolimod (orb2655858)

Catalog Number	orb2655858
Category	Small Molecules
Description	An immunomodulatory agent can be phosphorylated by SphK, then functions as a potent agonist of S1P receptors (S1P1/3/4/5) with EC50 of 0.3-5 nM.
Target	Lysophospholipid Receptor
Purity	>98% (HPLC)
MW	307.5
CAS Number	162359-55-9
Formula	$C_{19}H_{33}NO_2$
Storage	Storage temperature: -20°C. Stability: $\geq$ 2 years
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)