

Product Datasheet

Thalidomide-alkyne-C4-NHBoc (orb2663560)

Description	Thalidomide-alkyne-C4-NHBoc serves as an E3 Ligase Ligand-Linker Conjugate. It is utilized in the synthesis of PROTAC JAK1 degrader 1.
Storage	-20°C
Note	For research use only
Formula	$C_{24}H_{27}N_3O_6$
Purity	98.00%
MW	453.49
CAS Number	2093387-85-8
SMILES	<chem>O=C(OC(C)(C)C)NCCCCC#CC1=CC=CC=CC(=O)N(C(=O)C12)C3C(=O)NC(=O)CC3</chem>
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)