

Product Datasheet

Lenalidomide-acetylene-C3-MsO (orb2663906)

Description	Lenalidomide-acetylene-C3-MsO is a conjugate of an E3 ligase ligand and a linker (E3LigaseLigand-Linker Conjugates), utilized in the synthesis of SIAIS039.
Storage	-20°C
Note	For research use only
Formula	$C_{19}H_{20}N_2O_6S$
Purity	98.00%
MW	404.44
CAS Number	2138441-31-1
SMILES	<chem>O=C1C=2C=CC=C(C#CCCCOS(=O)(=O)C)C2CN1C3C(=O)NC(=O)CC3</chem>
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)