

Product Datasheet

L18I (orb2665228)

Description	L18I is a PROTAC (Proteolysis Targeting Chimera) that targets Btk to mitigate inflammation in autoimmune diseases. The components of L18I include the PROTAC target protein ligand IBT6A, the PROTAC linker Propargyl-PEG3-alcohol, and the E3 ubiquitin ligase ligand Lenalidomide-Br. Notably, the active control for the target protein ligand is IBT6A-CO-ethyne, while the conjugate of the E3 ubiquitin ligase ligand with the linker is Lenalidomide-C3-PEG3-N3.
Storage	-20°C
Note	For research use only
Formula	$C_{47}H_{51}N_{11}O_8$
Purity	98.00%
MW	897.98
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)