

Product Datasheet

Thalidomide-piperidine (orb2665541)

Description	Thalidomide-piperidine is an E3 ubiquitin ligase ligand-linker conjugate (E3 Ligase Ligand-Linker Conjugates) that can be utilized in the synthesis of SMART1.
Storage	-20°C
Note	For research use only
Formula	$C_{18}H_{19}N_3O_4$
Purity	98.00%
MW	341.36
CAS Number	2924858-38-6
SMILES	<chem>O=C1NC(=O)C(N2C(=O)C3=CC=CC(=C3C2=O)N4CCCCC4)CC1</chem>
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)