

## Product Datasheet

### Anti-CD11a/LFA-1 $\alpha$ Antibody (R7-1) (orb2665707)

<b>Description</b>	The Anti-CD11a/LFA-1 $\alpha$ Antibody (R7-1) is a chimeric mouse IgG1 antibody specifically targeting human CD11a/LFA-1 $\alpha$ . The recommended isotype control for this antibody is Mouse IgG1 kappa, Isotype Control.
<b>Conjugation</b>	Unconjugated
<b>Target</b>	Anti-CD11a/LFA-1 $\alpha$ Antibody (R7-1)
<b>Storage</b>	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
<b>Note</b>	For research use only
<b>Purity</b>	98.00%
<b>Expiration Date</b>	12 months from date of receipt.

---

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)