

## Product Datasheet

### Anti-GC1q R/C1QBP Antibody (60.11) (orb2665720)

|                        |                                                                                                                                                                                                               |
|------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Description</b>     | Anti-GC1q R/C1QBP Antibody (60.11) is a mouse IgG1, $\kappa$ chimeric antibody designed to target human GC1q R/C1QBP. The recommended isotype control for this antibody is Mouse IgG1 kappa, Isotype Control. |
| <b>Conjugation</b>     | Unconjugated                                                                                                                                                                                                  |
| <b>Target</b>          | Anti-GC1q R/C1QBP Antibody (60.11)                                                                                                                                                                            |
| <b>Storage</b>         | Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.                                                                       |
| <b>Note</b>            | For research use only                                                                                                                                                                                         |
| <b>Purity</b>          | 98.00%                                                                                                                                                                                                        |
| <b>Expiration Date</b> | 12 months from date of receipt.                                                                                                                                                                               |

---

#### Biorbyt Ltd.

7 Signet Court, Swann's Road,  
Cambridge, CB5 8LA, United Kingdom  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

#### Biorbyt LLC.

68 TW Alexander Drive,  
Durham, NC, 27713, United States  
Email: [info@biorbyt.com](mailto:info@biorbyt.com), [support@biorbyt.com](mailto:support@biorbyt.com)  
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)