

Product Datasheet

PE-conjugated Anti-CXCR2 antibody(DMC681); IgG1 Chimeric mAb (orb2669425)

Description	PE-conjugated Anti-CXCR2 antibody(DMC681); IgG1 Chimeric mAb
Species/Host	Rabbit
Reactivity	Human
Conjugation	Unconjugated
Tested Applications	FC
Target	CXCR2
Storage	Store at 2°C-8°C for 6 months
Note	For research use only
Isotype	Rabbit/Human Fc chimeric IgG1
Clonality	Monoclonal
Clone Number	DMC681
Uniprot ID	P25025, Q53PC4
Dilution Range	Flow Cyt 1:100
Expiration Date	12 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)