

Product Datasheet

Actb antibody (Biotin) (orb26948)

Description	Actb antibody (Biotin)
Species/Host	Rabbit
Reactivity	Mouse
Conjugation	Biotin
Tested Applications	ELISA
Immunogen	Recombinant Mouse Actin, cytoplasmic 1 protein (11-375AA)
Target	Actb
Preservatives	Preservative: 0.03% Proclin 300. Constituents: 50% Glycerol
Form/Appearance	Liquid
Storage	Upon receipt, store at -20°C or -80°C. Avoid repeated freeze.
Note	For research use only
Isotype	IgG
Clonality	Polyclonal
Uniprot ID	P60710
Dilution Range	ELISA
Expiration Date	12 months from date of receipt.