

Product Datasheet

Human IL-33 Protein, His Tag (orb2695786)

Description	Human IL-33 (C208S, C227S, C232S, C259S) Protein, His Tag (IL3-H52H6) is expressed from human 293 cells (HEK293). It contains AA His 109 - Thr 270 (Accession # O95760-1(C208S,C227S,C232S,C259S)).
Reactivity	Human
Endotoxins	1.0 EU per µg
Conjugation	Unconjugated
Target	IL-33
Preservatives	PBS, pH 7.4
Form/Appearance	Powder
Storage	-20°C to -70°C for 12 months in lyophilized state;. -70°C for 3 months under sterile conditions after reconstitution. For long term storage, the product should be stored at lyophilized state at -20°C or lower
Tag	C-10×His
Note	For research use only
Application notes	This protein carries a polyhistidine tag at the C-terminus. The protein has a calculated MW of 20.2 kDa. The protein migrates as 25-30 kDa when calibrated against Star Ribbon Pre-stained Protein Marker under reducing (R) condition (SDS-PAGE) due to glycosylation.
Purity	90%
MW	20.2 kDa
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44201223859353) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+14159065211) | Fax: [+1 \(415\) 651-8558](tel:+14156518558)