

Product Datasheet

Mouse TNFRSF1B Protein, His Tag (orb2699400)

Catalog Number	orb2699400
Category	Proteins
Description	Mouse TNFRSF1B Protein, His Tag
Target	TNFRSF1B
Reactivity	Mouse
Purity	>97%, determined by SDS-PAGE
Expression System	HEK293 Cells
Endotoxins	>97%, determined by SDS-PAGE
Source	This protein contains the extracellular domain of mouse TNFRSF1B (NP_035740.2) (Met 1-Gly 258) was fused with a polyhistidine tag at the C-terminus and expressed from HEK293 Cells.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)