

Product Datasheet

Human PROS1 Protein, His Tag (orb2700453)

Catalog Number	orb2700453
Category	Proteins
Description	Human PROS1 Protein, His Tag
Target	PROS1
Reactivity	Human
Purity	>95%, determined by SDS-PAGE
Expression System	HEK293 Cells
Endotoxins	>95%, determined by SDS-PAGE
Source	This protein contains the human PROS1 (P07225) (Met1-Ser676) was expressed with a polyhistidine tag at the C-terminus and expressed from HEK293 Cells.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)