

Product Datasheet

Human IL5RA Protein, His Tag (orb2701557)

Catalog Number	orb2701557
Category	Proteins
Description	Human IL5RA Protein, His Tag
Target	IL5RA
Reactivity	Human
Purity	>85%, determined by SDS-PAGE
Expression System	HEK293 Cells
Endotoxins	>85%, determined by SDS-PAGE
Source	This protein contains the human IL5Ra isoform 2 (Q01344-2) (Met 1-Arg 335) was fused with a polyhistidine tag at the C-terminus and expressed from HEK293 Cells.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com

Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)