

Product Datasheet

Human IL1R1 Protein, His Tag (orb2701864)

Description	Human IL1R1 Protein, His Tag
Reactivity	Human
Endotoxins	> 97%, determined by SDS-PAGE
Target	IL1R1
Note	For research use only
Purity	> 97%, determined by SDS-PAGE
Source	This protein contains the extracellular domain of human IL1R1 (NP_000868.1) (Met 1-Thr 332) was expressed with a fused polyhistidine tag at C-terminus and expressed from HEK293 Cells.
Expiration Date	6 months from date of receipt.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)