

Product Datasheet

Etokimab Biosimilar (Anti-IL-33 Reference Antibody) (orb2703879)

Catalog Number	orb2703879
Category	Antibodies
Description	Etokimab (Antibody ANB 020) is a humanized monoclonal antibody that targets IL-33. Etokimab can be used for the research of atopic dermatitis.
Target	Etokimab BiosimilarAnti-IL-33 Reference Antibody
Clonality	Monoclonal
Conjugation	Unconjugated
Antibody Type	Biosimilar Antibody
CAS Number	2022981-44-6
Storage	Maintain refrigerated at 2-8°C for up to 2 weeks. For long term storage store at -20°C in small aliquots to prevent freeze-thaw cycles.
Note	For research use only

Biorbyt Ltd.

7 Signet Court, Swann Road
Cambridge
CB5 8LA
United Kingdom

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\)1223 859353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)

Biorbyt LLC

68 TW Alexander Drive
Research Triangle Park
Durham
NC 27713
United States

Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)906-5211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)651-8558)